

Millennium B μ

Small dimensions... Big performance!

To offer the maximum in safety, speed, flexibility and versatility for your daily sterilization activities, Mocom has designed Millennium Bmicro, the small sized vacuum sterilizer (5.5 litre chamber) for rapid sterilization of all instrument loads in the dental practice.

All Millennium Bmicro sterilization cycles are equipped with pre-vacuum air removal to provide faster and more positive steam penetration throughout the entire load. The post-vacuum drying of wrapped and unwrapped instruments facilitates easy handling and sterile storage of instruments (wrapped only). Cycle selection, via the keypad, is simple and straight forward for the operator. An optional external printer records the performance of every cycle.



The Millennium Bmicro covers a wide range of situations and applications, adapting to the needs of all users. It is the fastest sterilizer from Start-to-Dry available in the market today:

only
18
minutes

Only 18 minutes to sterilize and dry wrapped hollow (handpieces) and solid loads!

only
9
minutes

Only 9 minutes to sterilize and dry unwrapped solid instruments in the Emergency cycle!

mocom
Sistemi innovativi di sterilizzazione

Quality System certified by:
CSQ CSQ MED